

# Invia® White Foam NPWT

## NEGATIVE PRESSURE WOUND THERAPY

- EN Instructions for use
- ES Instrucciones de uso
- FR Mode d'emploi
- DE Gebrauchsanweisung
- IT Istruzioni per l'uso
- SV Bruksanvisning
- DA Brugsanvisning
- FI Käyttöohjeet
- NO Bruksanvisning
- NL Gebruiksaanwijzing

## Table of contents

Intended use .....	4
Dosage form.....	4
Indications.....	4
Contraindications .....	4
Composition and characteristics .....	5
Operating instructions .....	5
Warnings .....	6
Precautions .....	6
Signs and symbols .....	7

## Intended use

The Invia White Foam Dressing is a polyvinyl alcohol (PVA) foam for Negative Pressure Wound Therapy (NPWT). The Invia White Foam Dressing is intended for use in conjunction with the Invia Negative Pressure Wound Therapy (NPWT) Systems.

## Dosage form

The Invia White Foam Dressing is supplied in a sterile packaging.

## Indications

The Invia White Foam Dressing can be used for almost all skin injuries to accelerate the formation of granulation tissue.

The Invia White Foam is used in combination with the Invia Negative Pressure Wound Therapy (NPWT) Systems.

Invia White Foam Dressing is used for:

- Traumatic wounds with delayed primary closure.
- Surgical wounds with delayed primary closure.
- Complex wounds (subverted wounds, tunnel wounds, wounds with exposed bone).
- Infected wounds with or without exposed bone, after debridement.
- Sternum wounds after median sternotomy without the risk of direct contact with the heart or larger vascular tissue.
- Superficial burns or scalds where a granulation into the foam is not desired.
- Fixation of free or pediculated lappets, split skin and full thickness skin grafts.
- Fixation of skin replacement material, skin allografts, artificial skin, cultivated and cultured skin tissue.
- Fistula, but only after surgical inspection and/or exploration.
  - Leg ulcers (ulcus cruris venosum) as a result of venous insufficiency.
- Diabetic foot ulcers, with or without neuropathy.
- Pressure ulcers, after surgical debridement of necrotic tissue.
- Pilonidal sinus.
- Pyoderma gangraenosum.
- Palliative wound care, but not for malignant wounds as a result of cancer or metastases.

## Contraindications

Do not use the Invia White Foam in the following cases:

- Necrotic wounds without previous debridement.
- Malignant wounds as a result of cancer or metastases.
- Unknown origin fistulas without previous surgical inspection and/or exploration.
- Increased bleeding risk of patients who receive anticoagulant agents or suffer from a coagulopathy as well as vascular anastomosis.
- Osteomyelitis.

## Composition and characteristics

The Invia White Foam is made from 100% polyvinyl alcohol.

The Invia White Foam is a hydrophilic, enormously absorbent foam. This means it is ideal for the effective removal of wound exudate and tissue fragments.

Combined with the Invia Negative Pressure Wound (NPWT) Therapy Systems with negative pressure levels between 80 mmHg and 150 mmHg, Invia White Foam ensures reliable removal of wound exudate from the wound. The applied negative pressure is distributed evenly on the base of the wound.

An intact sealing system guarantees a reliable airtight closure of the wound and thus an optimal use of negative pressure and exudate transport.

## Operating instructions

- Invia White Foam must only be used by medical and nursing staff.
- Take Invia White Foam from undamaged packaging.
- If the packaging is damaged, you must not use Invia White Foam.
- Follow all the usual guidelines for sterile treatment when applying Invia White Foam.
- Always cut Invia White Foam with sterile instruments to the desired size.
- Make sure that no scraps come into the wound bed when cutting Invia White Foam.
- Injuries of the deep tissue can be filled with several layers of Invia White Foam.
- If necessary, Invia White Foam can be used with appropriate antiseptics, either prior to fitting or using installation during the negative pressure therapy.
- Record the number of Invia White Foam pieces brought in the wound cavity in the patient's chart.
- Do not apply alcoholic solutions, as this compromises the integrity of Invia White Foam and its functioning.
- Make sure that Invia White Foam does not cover the skin around the wound or the growing epithelium.
- If there is evidence of an infection of the wound, Invia White Foam must be changed regularly, at least once every 24 hours until the signs of the infection have gone down significantly or disappeared altogether.
- An intermediate, non-stick wound cushion (limiting membrane) must be used when Invia White Foam comes in close contact with (parts of) blood vessels, tendons, nerves and other vulnerable vital organs and tissues.
- Always look for unexpected blood loss and instruct the patient accordingly.
- In the treatment of sternum wounds, controlling the blood pressure and ventilation is mandatory.
- Invia White Foam must be changed regularly within a maximum of 72 hours.
- If you see any signs of dehydration and adherence, Invia White Foam must be rehydrated with Ringer's solution first before it can be removed from the wound bed.
- Never use previously used Invia White Foam again.
- Discard Invia White Foam in accordance with local regulations.

## Warnings

Remove Invia White Foam in any case, if the patient shows unexpected symptoms such as:

- Increase in pain.
- Excessive accumulation of wound exudate.
- Clear signs of blood loss.
- Faecal or other excrements extravasation.
- In sternum wounds possibly decreased or impaired cardiac function.
- Clear signs of an infection.
- Suspicion of sepsis.
- Unexpected or unexplained deterioration of the general condition.

## Precautions

- Do not use Invia White Foam, if the packaging is damaged.
- Invia White Foam must not be re-sterilized.
- Do not dampen Invia White Foam with alcoholic solutions because this compromises the composition and integrity of the foam.
- Do not cut Invia White Foam above the wound, so that no foam fragments fall in the wound bed.
- If Invia White Foam is cut into pieces and placed in the wound, all pieces must be carefully counted and recorded in the patient's chart.
- When using Invia White Foam, a regular examination of the wound and also the unaltered condition of the patient is advisable.
- When using Invia White Foam in infected wounds, more frequent investigations of the wound and/or dressing changes are necessary.
- Invia White Foam is only used after proper wound debridement.
- After curing and adherence, Invia White Foam must be rehydrated and soaked first before it can be removed from the wound bed.

## Signs and symbols

This symbol indicates to follow instructions for use.

This symbol indicates the temperature limitation for operation transport and storage.

This symbol indicates the manufacturer.

This symbol indicates that the device should not be used after the date shown.

This symbol indicates to keep product dry.

This symbol indicates manufacturer's batch code.

This symbol indicates a single use device. Do not reuse the device.

This symbol indicates the device is in conformance with the Medical Device Directive 93/42/EEC (for EU only).

This symbol indicates the date of manufacture.

This symbol indicates the device is sterilized using irradiation.

This symbol indicates a carton package.

This symbol indicates to keep product away from sunlight.

This symbol indicates do not use the device if package is damaged.

This symbol indicates manufacturer's catalogue number.

This symbol indicates that the device should not be re-sterilized after once has been sterilized.

This symbol reminds consumer to dispose waste responsibly.

This symbol indicates number of items.

This symbol indicates that the material is part of a recovery/recycling process.

This symbol indicates that the material is part of a recovery/recycling process.

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Medela AG